



### **Full Cast Crown (Non-**

The following is about Full Cast Crown (Non-precious) related, I hope to help you better understand Full Cast Crown (Non-precious) .

### **Full Cast Crown (Non-precious ) Outsourcing Dental Lab in China,Manufacturer With Lab to Lab Service**



**Introduction:** FCC Nickel Free

**Material:** Argeloy N.P. Supreme Alloy

#### **Material List**

Type	Name of material	CE No.	Lot No.	Manufacturer	Composition
Non-precious Alloy	Argeloy N.P.Superme	CE0086	14034132	Argen(USA)	Co60%,Cr27%,W5%,Mo6%
Precious Alloy	Argedent 74	CE0086	718846622	Argen(USA)	Au74%,In2.5%,Pt4.0%,Pd7.6%
	Argedent 52	CE0086	71977557	Argen(USA)	Au52.5%,Pd26.9%,Ag16%,In2.5%
Semi-precious	Argistar 43+	CE0086	7545643	Argen(USA)	Au2%,Pd42.9%,Ag46%,In8%
Metal Framework	Argeloy N.P.Partial	CE0086	13869812	Argen(USA)	Co64%,Cr28%,Mo6%,
	Vitallium 2000	CE01663	N001830	Dentsply(USA)	Co63.8%,Cr28%,Mo5.9%
Acrylic Resin	TCS Unbreakable	CE0470	73021314	TCS(USA)	Acrylic Resin
	Vertex	CE120	XT016907	Castaperess	Acrylic Resin
	Lucitone 199	CE PR017600	131114	Dentsply(USA)	Acrylic Resin
All Ceramic	E max Press 3	CE0123	P19900	Ivoclar(Liechtenstein)	SiO279%,LiO15%,K2O5%,P2O5%,ZrO2%,ZrO3%
	Zirconia	HD601143840001	W1802260NG-1	Aidite(China)	ZrO2,Y2O3,AL2O3
	Full Contour Zirconia	HD601143840001	W171023QGC-1	Aidite(China)	ZrO2,Y2O3,AL2O3
Ceramic Powder	E max Ceram 3	CE0123	P62370	Ivoclar(Liechtenstein)	SiO265%,Al2O38%,Na2O6%,ZrO3%
	Vita Ceram	CE0124	LOT45190	Vita(Germany)	SiO265%,Al2O315%

We pursue the management tenet of "Quality is top-quality, Company is supreme, Status is first", and will sincerely create and share success with all shoppers for China New Design Full Cast Crown (Non-precious), Our Firm Core Principle: The prestige initially; The quality guarantee; The customer are supreme.

China New Design China Dental Lab, Good Quality and Good Esthetics, Strong infrastructure is the need of any organization. We are backed with a robust infrastructural facility that enables us to manufacture, store, quality check and dispatch our products worldwide. To maintain smooth work flow, we have sectioned our infrastructure into a number of departments. All these departments are functional with latest tools, modernized machines and equipment. Owing to which, we are able to accomplish voluminous production without compromising upon the quality.