




High Noble

The following is about PFM-Porcelain Fused to Metal High Noble related, I hope to help you better understand PFM-Porcelain Fused to Metal High Noble.

China High Noble outsourcing dental lab, manufacturer and Dental implants dental center with lab to lab service



Introduction : PFM-Porcelain Fused to Metal High Noble

Material: Argendent Euro (CE2797) /Vita vmk95

Composition of Metal(%)

Argendent Euro (CE2797) :Au 40% Pd 39.4% Ag 10% Ga
1.4% Ru Sn <1%

Argentent AH (CE0086) : Au 86.1% Pt 8.5% Pd 2.6% In
1.4% Ir Fe Ta<1%

